

Лазер	Channel	Фильтр	Красители
488	BL1	530/30	FITC, Alexa Fluor 488, CFSE, Fluo-3
	BL2	577/25	PE, PI
	BL3	620/30	ECD, PE-Texas Red, PE-CF594, PI
	BL4	690/50	PerCP, PC5, PC5.5, PerCP-Cy5.5, PI
	BL5	780/60	PE-Cy7
638	RL1	660/20	APC, Alexa Fluor 647, eFluor 660
	RL2	710/20	PC-Alexa Fluor 700, Alexa Fluor 700
	RL3	780/60	APC-Cy7, APC-Alexa Fluor 750, APC-H7, APC-eFluor 780
405	VL1	450/50	Pacific Blue, V450, eFluor 450, BV421
	VL2	530/30	Krome Orange, AmCyan, V500, BV510
	VL3	610/20	Violet610, BV605, Qdot 605
	VL4	660/20	Violet660, BV650, Qdot 655
	VL5	710/50	BV711, Qdot 705
	VL6	780/60	Violet780, BV785, Qdot800